



## In-situ gas analysis for fermentation

- or the easy way to control my process!

Tel. +49 - 2366 - 305-301 | Fax +49 - 2366 - 305-300 | Kontakt@BlueSens.de | www.BlueSens.de | Gas sensors „Made in Germany“



## BlueSens gas sensor GmbH

- Founded in July 2001 by Dr. Udo Schmale and Dr. Holger Mueller
  - up to 11 employees in 2009
  - global network of distributors
  - market leader for gas sensors in bioprocesses in Europe
- Main Focus: Gas analysis to monitor bioprocesses
  - Pharmaceutical industry
  - Research and development
  - Food industry
  - regenerative energy (biogas)
- Philosophy: Measuring where something happens (no sampling)



## Main Market: Pharmaceutical Industry

- Process optimization (R&D)
- controlled Scale up
- Quality control (yeast, yogurt etc.)
- Process control
  - Start culture (shake flask) to determine the optimal time for inoculation
  - Fermentation process (e.g. insulin production, cell culture, yeast production)



## Biogas Market

### Laboratory scale yield tests

- CH<sub>4</sub> quality
- Produced gas volume

### Biogas plant


- CH<sub>4</sub> control just in the near of the engine
- Total analysis of CO<sub>2</sub>, CH<sub>4</sub>, O<sub>2</sub> and H<sub>2</sub>S



## Laboratory Scale




## Gas analysis for the Scale up: “normal scale up”





## Gas analysis for the Scale up: "modern scale up"



- Online calculation of OUR, CER and RQ
- Increasing reproducibility
- Increasing productivity
- Greater understanding about your process
- Metabolic flux analysis and mass-balance calculations




## Gas analysis for the Scale up: "modern scale up"






## Principal of BCpreFerm



## Example 1

Optimisation of media with BCpreFerm




Source: Bayer Healthcare AG, Wuppertal





## Example 2





Source: Hexal Forschungs GmbH



## Example 3

Carbon balance in *S. cerevisiae* by online and offline analysis





## Example 3

### Physiological growth parameter and carbon balance

Spec. rates [mmol /gh]		Spec. yield coefficient [Cmol /Cmol]	
$r_{\text{glucose}}$	20.21 ± 0.61	$Y_{\text{glucose}}$	7.55 ± 0.16
$r_{\text{glycerol}}$	2.04 ± 0.34	$Y_{\text{glycerol}}$	0.38 ± 0.06
$r_{\text{acetate}}$	0.98 ± 0.09	$Y_{\text{acetate}}$	0.12 ± 0.01
$r_{\text{ethanol}}$	29.95 ± 0.66	$Y_{\text{ethanol}}$	3.73 ± 0.09
$r_{\text{CO}_2}$	34.37 ± 0.70	$Y_{\text{CO}_2}$	2.14 ± 0.03
		$\Sigma C$	-0.18 ± 0.06
		$\Sigma C/Y_{\text{gluc}}$ [%]	-2.36 ± 0.81

Carbon balance could be exactly closed and is reproducible

⇒ good method for quantitative physiological studies




## Production scale





## Control of fermentation processes

“normal” gas analysis




- a lot of maintenance
- delay of data
- difficult to control the process

Transport of the gas to the analyzer!



## Process control with in-situ gas sensors

“modern” gas analysis conform to PAT (process analytical technology)



- low maintenance
- online data
- easy to control the process

Measuring where something happens!





## Process control with in-situ gas sensors



for small  
pipes...



...for huge  
pipes...



...or for  
tubes

The size doesn't  
matter!




## Advantages of BlueSens' sensors

- small investment cost
- no sampling
- no gas treatment
- very stable and reliable measurements
- usable with every fermenter brand
- high redundancy through multiple sensors
- easy installation
- fits to every tube or pipe




## Process control with multiple sensors

- easy communication via Ethernet
- complete control loops
- different drivers available



## Example 4

Measurement of the Oxygen uptake rate (OUR) and Carbon dioxide emission rate (CER) of a yeast fermentation



Source: Schering AG, Engineo GmbH, BlueSens gas sensor GmbH



## Example 6

Discontinuous feeding of molasses based on an ethanol control

